



GetMiu V1.00

Data processing software for Hyperspectral Imager
powered by MATLAB R2017a

Features:

- Quasi fermi-level splitting (QFLS) mapping
- Urbach energy (Eu) mapping
- Voc prediction

Running environment

- Screen resolution: >1920*1080




Outline

- Running procedures (p3)
- Overview of the interface (p4)
- Data load (p5)
- Data selection (p6-p8)
- PL curve fit - 1D (p9-p21)
- PL curve fit - 2D (p22-26)
- Voc prediction (p27-p28)




Running procedures

Step 1: Copy the files below to a new folder



Step 2: open **GetMiu.m** with MATLAB, then you can see the raw codes in the Editor window (right figure→)


Step 3: Press **F5** or click the 'Run' button located in the menu bar to run the codes



```
Editor - GetMiu.m
fitdiff.m x len.m x Dalpha.m x TAlload.m x TAs.m x GetMiu.m x +
1 function varargout = GetMiu(varargin)
2 % GETMIU MATLAB code for GetMiu.fig
3 % GETMIU, by itself, creates a new GETMIU or raises the existing
4 % singleton*.
5 %
6 % H = GETMIU returns the handle to a new GETMIU or the handle to
7 % the existing singleton*.
8 %
9 % GETMIU('CALLBACK',hObject,eventData,handles,...) calls the local
10 % function named CALLBACK in GETMIU.M with the given input arguments.
11 %
12 % GETMIU('Property','Value',...) creates a new GETMIU or raises the
13 % existing singleton*. Starting from the left, property value pairs are
14 % applied to the GUI before GetMiu_OpeningFcn gets called. An
15 % unrecognized property name or invalid value makes property application
16 % stop. All inputs are passed to GetMiu_OpeningFcn via varargin.
17 %
18 % *See GUI Options on GUIDE's Tools menu. Choose "GUI allows only one
19 % instance to run (singleton)".
20 %
21 % See also: GUIDE, GUIDATA, GUIHANDLES
22
23 % Edit the above text to modify the response to help GetMiu
24
25 % Last Modified by GUIDE v2.5 06-May-2020 11:33:50
26
27 % Begin initialization code - DO NOT EDIT
28 gui_Singleton = 1;
29 gui_State = struct('gui_Name',       mfilename, ...
30                  'gui_Singleton',   gui_Singleton, ...
31                  'gui_OpeningFcn', @GetMiu_OpeningFcn, ...
32                  'gui_OutputFcn',  @GetMiu_OutputFcn, ...
33                  'gui_LayoutFcn',  [], ...
34                  'gui_Callback',    []);
35
36 if nargin && ischar(varargin{1})
37     gui_State.gui_Callback = str2func(varargin{1});
38 end
```



Overview of the interface



- A:** Read the PL-file and absorbance data
- B:** Select the 2D PL data (Displayed in G1)
- C:** Curve fit for Single point PL (Displayed in G2)
- D:** Curve fit for 2D area PL (Displayed in G3)
- E1:** View the fitted result
- E2:** View, analyze and copy the fitted result
- F:** Voc prediction (the value of PLQE is needed)




Data load

```
Editor - P:\work\code\Matlab codes\diff fit\Esma\GetMiu.m
ftdiff.m x len.m x Delpham x Tload.m x TAs.m x GetMiu.m x +
1 function varargout = GetMiu(varargin)
2 % GETMIU MATLAB code for GetMiu.fig
3 % GETMIU, by itself, creates a
4 % singleton*.
5 %
6 % H = GETMIU returns the handle
7 % the existing singleton*.
8 %
9 % GETMIU('CALLBACK',hObject,event
10 % function named CALLBACK in GET
11 %
12 % GETMIU('Property','Value',...)
13 % existing singleton*. Starting
14 % applied to the GUI before GetM
15 % unrecognized property name or
16 % stop. All inputs are passed t
17 %
18 % *See GUI Options on GUIDE's To
19 % instance to run (singleton)".
20 %
21 % See also: GUIDE, GUIDATA, GUIHANDLE
22
```


Load the PL file:

Step 1: Click 'PL data path' to select PL file (top figure)

Step 2: Choose the *.h5 file below and wait for the process to be completed (~10 seconds)

 S49_surf_532 nm_325mW_1_10sec(cal).h5

Step 3: Now you can see the 2D mapping and the center point of the PL data (right figure→)





Data selection

- For each dimension (wavelength, x-axis and y-axis), you can reset data range and data step

Data initialization

From To ☐ Average the step


Wavelength (nm) 650 850 Step (pix) 1

X (pix) 1 1392 X step (pix) 1


Y (pix) 1 1040 Y step (pix) 1

Colormap jet

- Edit range here
→ Function: initialize the raw data



- Default setting
- Center point located at X=696, Y=520



- Set wavelength: 660-590
- Set X: 1-1000
- Set Y: 1-1000
- Center point changes to X=500, Y=500



Data selection

- For each dimension (wavelength, x-axis and y-axis), you can reset data range and data step

Data initialization

From To ☐ Average the step

Wavelength (nm) 650 850

X (pix) 1 1392

Y (pix) 1 1040



Colormap jet

Step (pix) 1

X step (pix) 1

Y step (pix) 1

- Edit **data step** here
→ Function: remove pixel overlap,
reduce data size



- Set X-step: 20
- Set Y-step: 20




Data selection

- For each dimension (wavelength, x-axis and y-axis), you can reset data range and data step

Data initialization

Wavelength (nm)	From: 650	To: 850	<input type="checkbox"/> Average the step	Step (pix)	1
X (pix)	1	1392		X step (pix)	1
Y (pix)	1	1040		Y step (pix)	1
Colormap	jet				

- Choose your favorite **colormap** here
- If the data steps have been reset, you can **average all the stepped data** as the new data
→ Function: smooth the reduced data




- Tick 'Average the step'



PL curve fit -1D

- Select the **position** of the PL spectrum



Tick the checkbox then you can choose the position by mouse click


An alternative way is selecting the X-axis and Y-axis

Preview of the PL spectrum



PL curve fit -1D

- Select the **position** of the PL spectrum





When 'Using averaged PL' is ticked, the 1d PL spectrum is the averaged value of all the 2D PL spectra.

Preview of the PL spectrum in log-scale



PL curve fit -1D


- Select the **wavelength** of the PL spectrum





PL curve fit -1D

- Fit Options




Absolute PL intensity calibration

PL geometry: **Reflection** (default) or **Integrating sphere**

Reflection: CCD only collect a small solid angle of the PL photons from the top surface of the specimen

Collection cone





Integrating sphere: CCD collects all the PL photons by an integrating sphere (may not be the case of Hyperspectral imager)



PL curve fit -1D

- Fit Options



Numerical aperture of **objective lens**:

Objective lens: for PL detection (default value: 0.45, collection angle $\sim 27.1^\circ$)

How does the code do the flux calibration?

Assuming photons obey Lambertian distribution:

Correction factor for objective lens:

$$f_{lens} = 1 - \cos^2 \theta_{lens} \approx 0.2075$$


Calibrated PL Flux :

$$I_{PL} = I_{measured} / f_{lens}$$



PL curve fit -1D


Fit Options



Sample information

Select the absorption spectrum:

- Firstly you need to load the absorption file. Press this button to open [abs.dat](#)
- For your own absorption files please format it in *.dat, *.txt or *.mat



abs.dat - Notepad


File Edit Format View Help	
900	2487.564081
899	2498.342274
898	2514.05027
897	2521.917925
896	2531.765338
895	2544.588376
894	2552.491538
893	2562.383451
892	2574.272813
891	2585.189724
890	2592.146004
889	2601.100285
888	2607.076376

Wavelength value in cm^{-1}



PL curve fit -1D


Fit Options



Sample information

Select the absorption spectrum:


- Make sure the unit of absorption spectrum is cm^{-1}
- For a file with multi-columns, change the column number to switch the absorption spectra





PL curve fit -1D

- Fit Options



Sample information

Input the film thickness (default value: 640 nm)

Input the reflectivity at film's surface (reflectivity for excitation source, default value: 20%)




PL curve fit -1D

Fit Options

Blue curve: absorptivity calculated from absorption coefficient

Red curve: absorptivity extracted from the PL spectrum with the help of fitted Miu and T



Generalized Planck law (GPL): only works for PL's left branch (high energy side)

$$I_{PL}(\varepsilon) = \left\{ \frac{2\pi\varepsilon^2}{h^3c^2} \right\} \cdot \left\{ \frac{1}{\exp\left(\frac{\varepsilon - (E_F^{CB} - E_F^{VB})}{kT}\right) - 1} \right\} \cdot \{a(\varepsilon, E_F^{CB}, E_F^{VB}, T)\}$$

Where h is Planck's constant ($4.136 \times 10^{-15} \text{ eV}\cdot\text{s}$), c is the speed of light ($2.998 \times 10^8 \text{ m/s}$), k is the Boltzmann constant ($8.617 \times 10^{-5} \text{ eV/K}$), ε is the photon energy in eV , $a(\varepsilon, E_F^{CB}, E_F^{VB}, T)$ is the occupation and temperature dependent spectral absorptivity, and I_{PL} is the spectral photoluminescence emission flux in units of $\text{photons}/(\text{m}^2 \cdot \text{s} \cdot \text{eV})$.

The expression of spectral absorptivity:

$$a(\varepsilon) = \frac{(1 - R)(1 - \exp(-\alpha d))}{1 - R \exp(-\alpha d)}$$

Fitted result: QFLS (Miu) and Temperature (T)

Optional: fitting with fixed temperature

Optional: copy fitted curves to clipboard



Hybrid perovskite films approaching the radiative limit with over 90% photoluminescence quantum efficiency


Ian L. Braly^{1,4}, Dane W. deQuilettes^{2,4}, Luis M. Pazos-Outón³, Sven Burke^{2,7}, Mark E. Ziffer², David S. Ginger^{2,4} and Hugh W. Hillhouse^{1*}

Reducing non-radiative recombination in semiconducting materials is a prerequisite for achieving the highest performance in light-emitting and photovoltaic applications. Here, we characterize both external and internal photoluminescence quantum efficiency and quasi-Fermi-level splitting of surface-treated hybrid perovskite ($\text{CH}_3\text{NH}_3\text{PbI}_3$) thin films. With respect to the material bandgap, these passivated films exhibit the highest quasi-Fermi-level splitting measured to date, reaching $971 \pm 0.7\%$ of the radiative limit, approaching that of the highest performing GaAs solar cells. We confirm these values with independent measurements of internal photoluminescence quantum efficiency of $91.9 \pm 2.7\%$ under 1 Sun illumination intensity, setting a new benchmark for these materials. These results suggest hybrid perovskite solar cells are inherently capable of further increases in power conversion efficiency if surface passivation can be combined with optimized charge carrier selective interfaces.



PL curve fit -1D

- Fit Options



Blue curve: absorptivity calculated from absorption coefficient

Red curve: absorptivity extracted from the PL spectrum with the help of fitted Miu and T

GPL + DOS: only works for PL's left branch (high energy side)

$$I_{PL}(\varepsilon) = \left\{ \frac{2\pi\varepsilon^2}{h^3c^2} \right\} \cdot \left\{ \frac{1}{\exp\left(\frac{\varepsilon - (E_F^{CB} - E_F^{VB})}{kT}\right) - 1} \right\} \cdot \{a(\varepsilon, E_F^{CB}, E_F^{VB}, T)\}$$

Where h is Planck's constant ($4.136 \times 10^{-15} \text{ eV}\cdot\text{s}$), c is the speed of light ($2.998 \times 10^8 \text{ m/s}$), k is the Boltzmann constant ($8.617 \times 10^{-5} \text{ eV/K}$), ε is the photon energy in eV , $a(\varepsilon, E_F^{CB}, E_F^{VB}, T)$ is the occupation and temperature dependent spectral absorptivity, and I_{PL} is the spectral photoluminescence emission flux in units of $\text{photons}/(\text{m}^2 \cdot \text{s} \cdot \text{eV})$.

The expression of spectral absorptivity with DOS:

$$a(\varepsilon) = \frac{(1 - R)(1 - \exp(-\alpha d \cdot (f_v - f_b)))}{1 - R \exp(-\alpha d \cdot (f_v - f_b))}$$

Fitted result: QFLS (Miu) and Temperature (T)


Optional: fitting with fixed temperature

Optional: copy fitted curves to clipboard



PL curve fit -1D

- Fit Options



Optional: fit to the measured absorption coefficient to get E_g and E_u (Only available for the methods of GPL + Urbach, GPL + Urbach +DOS)

Urbach energy: the sub-band tail of absorption coefficient can be described by exponential decay:

$$\alpha(E) \propto \exp\left(\frac{E - E_g}{E_u}\right)$$

Note: E_g is an effective optical band gap (different from the energy gap)

Fitted result: absorption coefficient shows E_u is ~30 meV




PL curve fit -1D

- Fit Options

Blue curve: absorptivity calculated from absorption coefficient

Red curve: absorptivity extracted from the PL spectrum with the help of fitted M_{iu} , T , E_g and E_u



Generalized Planck law + Urbach: works for the full PL spectrum.

Left branch: fitted by GPL method

Right branch: fitted by Urbach method

Fitted result: PL spectrum shows E_u is ~ 13.85 meV




PL curve fit -1D

- Fit Options

Blue curve: absorptivity calculated from absorption coefficient

Red curve: absorptivity extracted from the PL spectrum with the help of fitted μ , T , E_g and E_u



Generalized Planck law + Urbach + DOS: works for the full PL spectrum.

Left branch: fitted by GPL method

Right branch: fitted by Urbach method

Fitted result: PL spectrum shows E_u is ~ 13.86 meV



PL curve fit -2D

- Fit Options

Second step: 2D Fitting

☒ Enable mouse click 2500 points

X1 (pix) 1 X2 (pix) 50

Y1 (pix) 1 Y2 (pix) 50

Fit miu_2D Abort

Number of fitting loops


To abort the fitting

Edit fitting area by mouse click

First click: set the top-left point

Second click: set the bottom-right point

The selected area is indicated by the dashed square
(available if mouse click is enabled)





PL curve fit -2D

- Fit Options

Second step: 2D Fitting

☒ Enable mouse click 2500 points

X1 (pix) X2 (pix)

Y1 (pix) Y2 (pix)

Fit miu_2D Abort

Press this button to start the 2D fit
The fitting method is in line with the 1D fit

Fitting progress

Fitting: 5.7%...

Data initialization

From To ☒ Average the step

Wavelength (nm) Step (pix)

X (pix) X step (pix)

Y (pix) Y step (pix)

Colormap

Second step: 2D Fitting

☒ Enable mouse click 2500 points

X1 (pix) X2 (pix)

Y1 (pix) Y2 (pix)

Fit miu_2D Abort

First step: 1D fitting

☐ Enable mouse click X (pix) Y (pix)

☒ Using area averaged PL ☒ Log scale of Y-axis

Fit from to copy fitted curves

PL geometry Numrical Aperture of lens

Absorption (/cm) Film thickness (nm) Refl.

Method Optional: Fit sub-band of Abs

Fit PL spectrum

Miu

T (K) ☐ fix

Eg (eV) Urbach (meV)

Fitted result

View miu_2D View T_2D

View Eg_2D View Urbach_2D

Optional: Get miu_device

View fitted absorptivity View fitted PL

Solar spectrum Absorbed flux: 0

PLQE_ext T (K)

Miu_SQ Err


Miu_rad Err

Miu_dev Err

Calculate




PL curve fit -2D



The fitted result shows $\mu \sim 1.225 \pm 0.015$ eV

- This result is more accurate than the one got in page 21
- The width of μ is related to film's spatial inhomogeneity




View the fitted result: here μ_{2D} is displayed

Data manipulation: here μ_{2D} is selected and displayed in histogram




PL curve fit -2D



The fitted result shows
 $E_u \sim 14.1 \pm 0.3$ meV

- This result is more accurate than the one got in page 21
- The width of E_u is related to film's spatial inhomogeneity



View the fitted result: here E_u _2D is displayed

Data manipulation: here E_u _2D is selected and displayed in histogram



PL curve fit -2D

A reference based on the study of 2D PL spectra

Energy &
Environmental
Science



COMMUNICATION

View Article Online
View Journal | View Issue



Cite this: *Energy Environ. Sci.*,
2016, 9, 2286

Received 14th February 2016,
Accepted 25th May 2016

DOI: 10.1039/c6ee00462h

www.rsc.org/ees

Quantification of spatial inhomogeneity in perovskite solar cells by hyperspectral luminescence imaging†

Gilbert El-Hajje,^{ab} Cristina Momblona,^c Lidón Gil-Escrig,^c Jorge Ávila,^c
Thomas Guillemot,^{de} Jean-François Guillemoles,^{b,d} Michele Sessolo,^c
Henk J. Bolink*^c and Laurent Lombez*^b

QFLS mapping





Fig. 3 (a) Absolute mapping of the quasi-Fermi level splitting derived from photoluminescence (a and b, top) and electroluminescence (c and d, bottom), for perovskite solar cells using PCBM (left, device A) or C₆₀ (right, device B) as the ETL. The Gaussian distribution of the signal has been overlaid with the color scale as a guide to the eye.



Voc prediction

Red curve: absorbed solar spectrum under AM G1.5

AM G1.5



Simulated PL spectra with different type of QFLS

miu_SQ: PL at Shockley-Queisser limit
miu_rad: PL at radiative limit
miu: measured PL

Parameters for Miu_dev caculation

PLQE is also termed as PLQY

$$\mu_{dev} = \mu_{rad} - k_B T |\ln PLQE_{ext}|$$

Miu_dev is highly dependent on PLQE and Temperature



Voc prediction

A reference based on QFLS fitting:
Use high external PLQE to explain high Voc

There are amendments to this paper

LETTER

doi:10.1038/nature25989

Maximizing and stabilizing luminescence from halide perovskites with potassium passivation

Mojtaba Abdi-Jalebi¹, Zahra Andaji-Garmaroudi¹, Stefania Cacovich², Camille Stavrakas¹, Bertrand Philippe³, Johannes M. Richter¹, Mejd Alsari¹, Edward P. Booker¹, Eline M. Hutter⁴, Andrew J. Pearson¹, Samuele Lilliu^{5,6}, Tom J. Savenije⁴, Håkan Rensmo³, Giorgio Divitini², Caterina Ducati², Richard H. Friend¹ & Samuel D. Stranks¹



Hybrid perovskite films approaching the radiative limit with over 90% photoluminescence quantum efficiency

Ian L. Braly^{1,4}, Dane W. deQuilettes^{2,4}, Luis M. Pazos-Outón³, Sven Burke², Mark E. Ziffer², David S. Ginger^{2*} and Hugh W. Hillhouse^{2,1*}

Reducing non-radiative recombination in semiconducting materials is a prerequisite for achieving the highest performance in light-emitting and photovoltaic applications. Here, we characterize both external and internal photoluminescence quantum efficiency and quasi-Fermi-level splitting of surface-treated hybrid perovskite (CH₃NH₃PbI₃) thin films. With respect to the material bandgap, these passivated films exhibit the highest quasi-Fermi-level splitting measured to date, reaching $97.1 \pm 0.7\%$ of the radiative limit, approaching that of the highest performing GaAs solar cells. We confirm these values with independent measurements of internal photoluminescence quantum efficiency of $91.9 \pm 2.7\%$ under 1 Sun illumination intensity, setting a new benchmark for these materials. These results suggest hybrid perovskite solar cells are inherently capable of further increases in power conversion efficiency if surface passivation can be combined with optimized charge carrier selective interfaces.